



DMS/DMS Intermediate

The dynamic, angle-adjustable barrel plate system






Table of Contents

	Pages
DMS Product Features	4-7
Surgical Techniques	
■ DMS Surgical Technique	8-17
■ DMS Intermediate Surgical Technique	18-29
DMS Product Range	
■ Implants	30-31
■ Instruments	32-35
■ Storage Trays	36-37
DMS Intermediate Product Range	
■ Implants	38-39
■ Instruments	40-43
■ Storage Trays	44-45

DMS

Dynamic, angle-adjustable barrel plate system

DMS

Coxal-femoral fractures

First acceptable results following femoral neck and pertrochanteric femur osteosynthesis were reported by Smith-Petersen (1931), using a 3-lamellae pin with three-star cross section. Further developments included the implants according to Thornton (1937), Jewett (1941) and McLaughlin (1947), which are rigidly connected to a plate applied laterally to the femur. Systematic further developments were the angle plates designed by Schneider and Müller (1957).

E. Pohl (1951) can claim credit for having developed the first non-interlocking connection between an intra-medullary force carrier and a lateral anchoring plate. The principle consists of a dynamic connection between a lag screw and a femoral “barrel

plate” that allows self-compression (sliding barrel principle).

This brilliant idea made it possible to avoid many of the previous complications (head perforation, pseudarthrosis, secondary dislocation). Modified implants were introduced in the United States by Callender (1967) as “Sliding Hip Screws” and by Richards (1971) as “Compression Screws”. An optimized follow-up model was developed by the AO (1979) on the basis of Pohl’s system. While maintaining the „sliding barrel principle“, additional rotational stability is ensured by form-fit (i.e. by using a hexagonal screw instead of a round one and flattening the barrel on two sides to provide corresponding sliding surfaces).



DMS Intermediate

Also based on Pohl's system, we have now newly developed a dynamic, infinitely adjustable, variable barrel plate for stabilizing proximal femoral fractures. If required, the KLS Martin implant (DMS) allows valgisation of the femoral head fragment after fixation of the fracture and before compressing the pertrochanteric fracture surfaces, particularly in the case of unstable fractures.

The greatest advantage of this implant lies in the fact that the required valgisation angle can be adjusted to the individual patient's needs by means of a worm gear.

The DMS Intermediate is a third smaller than the DMS standard plate and therefore better suited for the anatomy of children and short adults. At the same time, the advantages of intraoperative adjustment of the angle of the implant are fully retained.

K. K. DITTEL

Feature, function and benefit



DMS



DMS Intermediate




The Dynamic Martin Screw (DMS) was specially designed to stabilize proximal femoral fractures. The dynamic, angle-adjustable barrel plate offers perfect shaft congruence for patient-specific adaptation. The “sliding barrel” principle of the implant enables full weight bearing at an early stage. And thanks to the universal applicability of the plate, stockkeeping needs and corresponding costs are significantly reduced.

Alternatively to the DMS standard plate there is available the DMS Intermediate which is well suited for the anatomy of the child and small adults as well.

DMS – simple and dynamic!

Plate	Features	Benefit
	<ul style="list-style-type: none"> ■ Dynamic ■ Infinitely adjustable ■ Angle-adjustable 	<ul style="list-style-type: none"> ■ Ideal shaft congruence ■ Patient-specific adaptation → low complication rates → increased safety ■ Intraoperative flexibility ■ Secondary varus or valgus correction possible
	<ul style="list-style-type: none"> ■ Can be inserted with maximum precision ■ Optimal adaptation of the plate to the bone 	<ul style="list-style-type: none"> ■ Good biomechanical conditions ■ Excellent and reliable stability under load even in difficult cases ■ Compression and reduction after plate fixation
	<ul style="list-style-type: none"> ■ “Sliding barrel” principle 	<ul style="list-style-type: none"> ■ Fast restoration of full weight-bearing capacity
	<ul style="list-style-type: none"> ■ One barrel plate replaces up to 6 fixed-angle plates ■ Universal application 	<ul style="list-style-type: none"> ■ Cuts stockkeeping costs significantly

Lag screw

	<ul style="list-style-type: none"> ■ DMS Intermediate: forward and backward cutting lag screw 	<ul style="list-style-type: none"> ■ Easy implantation and safe explantation
---	--	---

Step by step
to optimal fixation

DMS indications



Medial femoral neck fractures
DMS 2-hole plate



Pertrochanteric femoral fractures
DMS 4-hole plate



Intertrochanteric reversed fractures
DMS 4-hole plate



Subtrochanteric fractures
DMS 4-hole plate


DMC indications



Condylar femoral fractures



Supracondylar femoral fractures





Surgical technique

DMS surgical technique

Pertrochanteric femoral fracture

Pages 10-17






Preoperative X-ray




Positioning the patient

The patient is placed on a radiolucent operating table in dorsal decubitus position.



1. Lateral surgical approach acc. to Bauer

Starting two finger widths proximally to the tip of the trochanter major, a straight, lateral skin incision is made, length 15 cm.




2. Countersinking the cortex of the femur

The 4.5-mm drill bit is used to countersink the femoral cortex.




4.5-mm drill bit



3. Positioning the guide wire

The guide wire is positioned using the aiming device, which can be adjusted between 135° and 150°.



4. Inserting the guide wire

The guide wire is then inserted under image intensifier control, making sure that it lies centrally in the mid-axis of the femoral head.

Notice:

Alternatively, the guide wire can be inserted free-hand using an appropriate soft-tissue sleeve.



Guide wire




Aiming device



Guide wire




Soft tissue sleeve



5. Image intensifier control

Correct placement of the guide wire is checked with the image intensifier.




6. Determining guide wire length

Once the guide wire has been correctly positioned, its length can be conveniently read off the scale of the measuring sleeve.




Measuring sleeve



7. Drilling the hole

After setting it to the measured value (-10 mm), the DMS combo reamer is drilled into the bone along the guide wire under image intensifier control until the cone of the third stage has fully entered the lateral cortex.



8. Tapping

Optionally using the centering sleeve and the T-handle, the tap is now screwed in to a point 10 mm away from the cortex.

The depth of the thread can be directly read off the mark on the centering sleeve.



Combo reamer


Step 1:



Centering sleeve


Tap

Step 2:



T-handle

Tap with centering sleeve





9. Inserting the lag screw

The length of the lag screw is identical with the set drilling depth.

To insert the lag screw, it is first attached to the screwdriver and the connector before it is screwed in with the safety inserter, the 11-mm centering sleeve and the T-handle.

Notice: In hard bone, the lag screw is inserted up to the first mark, in osteoporotic bone up to the last mark.




10. Inserting the plate

Once the lag screw has been positioned correctly, the handle with the safety inserter and the centering sleeve can be removed.

Now a plate of correct length can be passed over the screwdriver onto the lag screw.

Notice:

If the plate is not parallel to the longitudinal axis of the femur, the T-handle can be applied again to fine-adjust the screw by turning it clockwise.




11. Adjusting the plate

Once the plate is in the correct position relative to the femoral axis, it is adjusted with the worm gear as required for proper valgus or varus correction.

Using the screwdriver, the worm gear is operated until the plate attaches perfectly to the femur.

To ensure a secure seat, the plate impactor is used to precision-adjust the DMS plate on the femur.



12. Fixing the plate

To fix the DMS plate in position on the femur, 4.5-mm cortical screws are used.

For the plate hole located directly underneath the worm gear, a 6.5-mm cancellous screw can also be used for fixation of the lesser trochanter.




Screwdriver



Plate impactor



Screwdriver






Drill bit for 4.5-mm
cortical screws



Depth gauge



Drill bit guide



13. Inserting the compression screw

Postoperative X-ray check

In a last step, the fracture is compressed by inserting the DMS compression screw. In osteoporotic bone, compression paths of up to 6 mm can occur.

Following compression, the compression screw is removed.

Left picture: Fracture without compression

Right picture: Compressed fracture



Screwdriver

DMS Intermediate – step by step to optimal fixation

DMS Intermediate indications

Proximal femoral fractures in short adults



Medial femoral neck fracture
DMS Intermediate 2-hole plate



Pertrochanteric femoral fracture
DMS Intermediate 4-hole plate



Intertrochanteric reversed fracture
DMS Intermediate 4-hole plate



Subtrochanteric fracture
DMS Intermediate 4-hole plate

Displacement osteotomy for coxa valga correction



preoperative



postoperative, correction approx. 10°
Fixation with DMS Intermediate


Displacement osteotomy for coxa vara correction



preoperative



postoperative, correction approx. 20°
Fixation with DMS Intermediate



Surgical technique

DMS Intermediate surgical technique

Intertrochanteric varisation osteotomy

with a correction angle of 20°

Prof. Meiß

Pages 20-29






Preoperative planning

As part of the preoperative planning, the scope of the varus correction needs to be determined with precision.

On this basis, the DMS Intermediate plate with the correct angulation is then selected.

Detailed planning requires a pelvic X-ray view or an AP X-ray of the hip.




Preoperative planning


A paper tracing is made from the X-ray image. The femur is then cut out and the desired correction set after dissection in the osteotomy region.

This can be done with or without wedge removal (i.e. lateral opening or medial closing).

In addition, the prospective position of the DMS Intermediate plate before and after “osteotomy” is sketched in as well, using an X-ray template.



X-ray template




Positioning the patient

The patient is placed in dorsal decubitus position on a radiolucent operating table. Both legs are covered so that they are easily movable.

Besides, it is advisable to place a radiolucent cushion or folded towel underneath the sacrum in order to lift the pelvis and thighs a little.

If pelvic osteotomy is additionally planned for the same session, the iliac crest is covered as well up to costal arch level.




1. Approach – opening

Surgical access is made by means of a lateral longitudinal thigh incision from a point slightly distally of the tip of the trochanter to the proximal shaft of the femur.

The fascia is split along the course of the fibers, with the leg kept in a slightly abducted position to reduce fascia tension.

Thereafter, the trochanteric bursa is opened and the tissue retracted ventrally and dorsally.




2. Approach – opening

The intertrochanteric femoral region must now be exposed.

Notice:

To ensure good orientation, subperiosteal exposure of the trochanter is recommended until the proximal end of the femoral neck is ventrally palpable or visible.



3. Placing the guide wire

To mark the femoral neck antetorsion, a 2-mm-thick Kirschner wire is inserted into the femoral neck as centrally as possible under image intensifier control.

Notice:

It is important to maintain the predetermined angle to the shaft. The entry point of the guide wire should be located approx. 1 cm distally of the trochanteric apophyseal growth plate.



Guide wire




4. Placing the guide wire

The angle to the shaft must be set so that the intended varus correction can be achieved with the selected DMS Intermediate plate by reducing the angulation.

The guide wire is not pushed up to the subchondral lamella of the femoral head, but to a point just short of the epiphyseal plate of the femoral head.

Notice:

Image intensifier control in two planes is essential!



5. Setting the drill

Once positioned correctly, the length of the guide wire inside the bone is determined with the measuring sleeve.

The (3-stage) combo reamer is then set to the measured value minus 5 mm.

Notice:

When loosening the knurled nut, please note that this is a left-hand thread (to be turned clockwise!) and that the nut must be removed completely.



Measuring sleeve



Combo reamer




6. Pre-drilling

Pre-drilling is done via the guide wire under irrigation until the cone of the third stage is fully immersed in the lateral cortex.

Notice:

In any case, it is important to do the pre-drilling under image intensifier control to prevent the guide wire from being inadvertently pushed forward into the growth plate.



7. Tapping


The tap is inserted to the true drilling depth. This can be checked via the scale marks at the level of the lateral cortex.

Notice:

When screwing in the tap, be sure not to overturn it, as stripping (destroying) the thread inside the bone would jeopardize proper lag screw anchorage.



Combo reamer



T-handle


Tap



8. Inserting the lag screw

The length of the selected lag screw must match the drilling depth set on the combo reamer.

Using the lag screw connector, the screw is first fitted with the lag screw inserter sleeve before it is finally screwed in with the safety inserter for lag screw inserters and the T-handle.



9. Trial insertion of the plate

Now the appropriate DMS Intermediate plate can be inserted for testing purposes. This is done by pushing the cylindrical (barrel) part of the plate over the lag screw inserter onto the lag screw and attaching the plate to the shaft of the femur.

Notice:


Due to the hexagonal design of the lag screw, the plate may initially be out of line with the femoral shaft. In this case, attach the T-handle with the safety inserter again and adjust the lag screw by rotating it a little forward or backward.



10. *Adjusting the plate*

The plate must now be adapted accurately to the femoral shaft. This is done by adjusting its angulation via the adjusting screw, using the 2.5-mm Allen screwdriver.

To perform the osteotomy, the barrel plate is removed again while leaving the angular position unchanged for the time being. The lag screw remains in its final position.



11. *Osteotomy*


To mark the rotation, or plan the rotational correction, two Kirschner wires with a thickness of approx. 2.5 mm are inserted into the femur ventrally or ventrolaterally, proximally and distally of the planned osteotomy.



Screwdriver



Guide wire




12. Osteotomy

Starting just distally of the lag screw insertion hole, the osteotomy is performed towards medial in such a way that the largest part of the lesser trochanter is left in place.

Preplacement of a Kirschner wire that can be used as a guide for the distal saw cut has proved useful.

Completion of the saw cut is followed by a marked widening of the osteotomy gap using a bone spreader in order to release the soft tissues.

Then the fragments are put in position according to the planned varus correction.



13. Inserting the plate

Now the plate can be inserted. This may require slight correction with the plate impactor.

Thereafter, the compression screw is inserted and slightly tightened.


The angulation is reduced as appropriate by rotating the adjusting screw clockwise (1 turn = 7.5°).



Plate impactor



Screwdriver



14. Fixing the plate

The plate is finally fixed in place with self-tapping 3.5-mm cortical screws.

If possible, the first screw should be inserted into the proximal sliding hole in compression position so the fragments are impacted a little.

When performing an opening osteotomy, it is recommended to fill the lateral gap by using a spongiosaplasty.

The compression screw is usually removed again to achieve optimal dynamization.

Finally, the result is documented in two planes by image intensifier.

15. Wound closure

Thorough irrigation of the wound.

Good refixation of the vastus lateralis muscle if it has been detached by an L-shaped incision.

Then layer-by-layer wound closure with loose suturing of the trochanteric bursa. If necessary, use of 1 or 2 drains. Application of a long-leg (pelvis-leg-foot) spica cast is required only in exceptional cases, because correct implantation of the DMS Intermediate plate provides solid fixation allowing early range-of-motion exercises if not partial weight bearing.



Screwdriver

Drill bit for
3.5-mm cortical
screws

Depth gauge

Drill bit guide



Follow-up-treatment

On the 2nd or 3rd postoperative day, radiographic assessment:
AP pelvic view and axial view of the hip (Lauenstein position).

Begin with mobilization on about the 3rd postoperative day with range-of-motion exercises including use of a continuous passive motion device and transfer in a wheelchair in recumbent position. Alternatively, mobilization with a walking frame or forearm crutches (as a rule with partial weight bearing).

Discharge on the 5th to 7th postoperative day.

After six weeks, radiographic follow-up and decision about full weight bearing.

Implant removal after three to four months.

DMS product range
Implants

DMS plate

Angulation
85°-145°

34 mm
Barrel length

St 1
unit(s)



Holes	Item No.
1	25-181-01-05
2	25-181-02-05
4	25-181-04-05
5	25-181-05-05
6	25-181-06-05
8	25-181-08-05
10	25-181-10-05
12	25-181-12-05
14	25-181-14-05
16	25-181-16-05
18	25-181-18-05

DMS plate

Angulation
85°-145°

24 mm
Barrel length

St 1
unit(s)



Holes	Item No.
4	25-182-04-05
6	25-182-06-05
8	25-182-08-05

DMC plate

Angulation
95°-115°




24 mm
Barrel length

St 1
unit(s)



Holes	Item No.
2	25-183-02-05
4	25-183-04-05
6	25-183-06-05
8	25-183-08-05
10	25-183-10-05
12	25-183-12-05
14	25-183-14-05
16	25-183-16-05

Icon explanations

-  Steel
-  Packaging unit
-  Hexagon head

Lag screws Ø 12.7 mm



Lag screws	Item No.
50 mm	25-180-50-05
55 mm	25-180-55-05
60 mm	25-180-60-05
65 mm	25-180-65-05
70 mm	25-180-70-05
75 mm	25-180-75-05
80 mm	25-180-80-05
85 mm	25-180-85-05
90 mm	25-180-90-05
95 mm	25-180-95-05
100 mm	25-180-00-05
105 mm	25-180-05-05
110 mm	25-180-10-05
115 mm	25-180-15-05
120 mm	25-180-20-05

Cortical screws Ø 4.5 mm



Lag screws	Item No.
30 mm	25-114-30-05
32 mm	25-114-32-05
34 mm	25-114-34-05
36 mm	25-114-36-05
38 mm	25-114-38-05
40 mm	25-114-40-05
42 mm	25-114-42-05
44 mm	25-114-44-05
46 mm	25-114-46-05
48 mm	25-114-48-05
50 mm	25-114-50-05
52 mm	25-114-52-05
56 mm	25-114-56-05
60 mm	25-114-60-05
64 mm	25-114-64-05
70 mm	25-114-70-05

Compression screw



25-180-99-05
Ø 4.5 mm
Lag screws 37 mm



DMS product range
Instruments



1/2

25-238-00-05
Guide wire
22.5 cm / 8 7/8"
Ø 2.5 mm
St 5 units



1/2

25-238-05-07
Measuring sleeve
20.5 cm / 8 1/8"
St 1 unit



1/2

25-238-13-07
Combo reamer
for lag screw and 34-mm barrel
22 cm / 8 5/8"
St 1 unit

25-238-12-07
Combo reamer
for lag screw and 24-mm barrel
St 1 unit






1/2

25-238-18-04
Centering sleeve
(11 mm) for 25-238-21-07
11 cm / 4 3/8"
Al 1 unit

25-238-19-04
Centering sleeve
for 25-238-06-07
Al 1 unit

Icon explanations

-  St Steel
-  Al Aluminum
-  1 unit(s) Packaging unit



25-238-17-07

Tap

23 cm / 9"



25-238-06-07

Lag screw inserter

21.5 cm / 8 4/8"



25-238-20-07

Lag screw connector

20 cm / 7 7/8"



25-238-21-07

*Safety inserter
for screwdriver*

24 cm / 9 4/8"

DMS product range Instruments

3.5 mm



25-238-25-04
Spare part, alone



1/2

25-268-35-07
Screwdriver

25 cm / 9 7/8"

for worm gears and for 4.5-mm
cortical screws and 6.5-mm
cancellous screws



1/2

25-238-14-07
Plate impactor

20 cm / 7 7/8"



1/2

25-238-08-07
T-handle

6.5 cm / 2 4/8"




1/2

25-212-45-07
Soft tissue sleeve

14.5 cm / 5 6/8"





90-651-51-21
X-ray template for
preoperative planning

1
unit(s)


Icon explanations

- St Steel
- 1 unit(s) Packaging unit
- Hexagon head



25-238-01-07
Aiming device

St 1
unit(s)



25-210-30-07
Drill bit
Ø 3.2 mm
with AO attachment
for 4.5-mm cortical screws


St 1
unit(s)



25-219-05-07
Depth gauge

St 1
unit(s)

For Ø 3.2-mm drill bits




25-222-44-01
Drill bit guide
16.5 cm / 6 1/2"
for 3.2-mm drill bits

St 1
unit(s)

DMS product range

Storage trays and set list

The instrument set described as follows represents the basic equipment needed for implanting DMS plates.



Set including options

KLS Martin DMS Set		
consisting of:		
Item No.	Quantity	Description
25-400-00-04		KLS Martin DMS Set
Implants		
25-180-60-05	1	Lag screw, 60 mm
25-180-70-05	1	Lag screw, 70 mm
25-180-80-05	1	Lag screw, 80 mm
25-180-85-05	1	Lag screw, 85 mm
25-180-90-05	1	Lag screw, 90 mm
25-180-95-05	1	Lag screw, 95 mm
25-180-00-05	1	Lag screw, 100 mm
25-180-10-05	1	Lag screw, 110 mm
25-180-99-05	3	Compression screws
25-181-02-05	1	Plate, 2-hole
25-181-04-05	1	Plate, 4-hole
25-181-05-05	1	Plate, 5-hole
25-181-06-05	1	Plate, 6-hole
Instruments		
25-238-00-05	1	Guide wire, 5 pcs.
25-238-05-07	1	Measuring sleeve
25-238-06-07	1	Lag screw inserter
25-238-08-07	1	T-handle
25-238-13-07	1	Combo reamer
25-238-14-07	1	Plate impactor
25-238-17-07	1	Tap
25-238-18-04	1	Centering sleeve (11 mm) 25-238-21-07
25-238-20-07	1	Lag screw connector
25-238-21-07	1	Safety inserter for lag screw inserter
22-368-35-07	1	Screwdriver, 3.5 mm, Allen
Storage		
55-550-60-01	1	Perforated tray for implants
55-555-61-04	1	Perforated tray for instruments
55-443-22-04	2	Coding labels, lettered, w/o hole
55-443-12-04	2	Logistics framelets, red
55-442-13-04	1	Container, 600 x 300 x 140 mm

DMS Intermediate product range
Implants

*DMS
Intermediate
plate*

Angulation
80°-110°

19 mm
Barrel length



STERILE !



Holes	Item No.
3	25-178-03-05
4	25-178-04-05
6	25-178-06-05
8	25-178-08-05

*DMS
Intermediate
plate*

Angulation
110°-145°

19 mm
Barrel length



STERILE !



Holes	Item No.
3	25-179-03-05
4	25-179-04-05
6	25-179-06-05
8	25-179-08-05



Icon explanations

St

Steel

1
unit(s)

Packaging unit

Hexagon head

Hexagon head

STERILE |

steam-sterilized

Lag screws

Ø 9.5 mm

St 1
unit(s)

STERILE |




Lag screws	Item No.
35 mm	25-170-35-05
40 mm	25-170-40-05
45 mm	25-170-45-05
50 mm	25-170-50-05
55 mm	25-170-55-05
60 mm	25-170-60-05
65 mm	25-170-65-05
70 mm	25-170-70-05
75 mm	25-170-75-05
80 mm	25-170-80-05
85 mm	25-170-85-05
90 mm	25-170-90-05

Cortical screws

Ø 3.5 mm

St 5
unit(s)

STERILE |



Lag screws	Item No.
16 mm	25-103-16-05
18 mm	25-103-18-05
20 mm	25-103-20-05
22 mm	25-103-22-05
24 mm	25-103-24-05
26 mm	25-103-26-05
28 mm	25-103-28-05
30 mm	25-103-30-05
32 mm	25-103-32-05
34 mm	25-103-34-05
36 mm	25-103-36-05
38 mm	25-103-38-05
40 mm	25-103-40-05
45 mm	25-103-45-05
50 mm	25-103-50-05
55 mm	25-103-55-05
60 mm	25-103-60-05

Compression screw



25-170-99-05


Ø 3.5 mm

Lag screws 25 mm

St 1
unit(s)

STERILE |

DMS Intermediate product range Instruments



25-210-25-07
Drill bit
Ø 2.5 mm
for 3.5-mm
cortical screws

St 1
unit(s)



25-219-04-07
Depth gauge
16 cm / 6 2/8"

St 1
unit(s)



25-222-32-01
*Drill bit guide for neutral
and eccentric holes*
14 cm / 5 4/8"

St 1
unit(s)



25-228-00-05
*DMS Intermediate
guide wire*
18 cm / 7"
Ø 2.0 mm

St 1
unit(s)

Icon explanations

St Steel

1 unit(s) Packaging unit



25-228-05-07
*DMS Intermediate
measuring sleeve*

15 cm / 6"

St 1 unit(s)



25-228-06-07
*DMS Intermediate
lag screw inserter*

7 cm / 2 5/8"

St 1 unit(s)



25-228-20-07
*DMS Intermediate
lag screw connector*

9 cm / 3 4/8"

St 1 unit(s)



25-228-21-07
*DMS Intermediate safety
inserter for lag screw inserter*

12.5 cm / 5"

St 1 unit(s)

DMS Intermediate product range
Instruments



25-228-12-07
*DMS Intermediate
combo reamer*
16 cm / 6 2/8"
St 1 unit(s)






25-228-14-07
*DMS Intermediate
plate impactor*
15 cm / 6"
St 1 unit(s)



25-228-17-07
*DMS Intermediate
tap*
15 cm / 6"
St 1 unit(s)

Icon explanations


-  Steel
-  Packaging unit
-  Hexagon head



25-238-08-07
DMS T-handle
12 cm / 4 7/8"



22-368-24-07
Screwdriver, Allen
24 cm / 9 4/8"
for worm gears and for cortical
screws Ø 3.5 mm and cancellous
screws Ø 4.0 mm



25-238-01-07
X-ray template for
preoperative planning





DMS Intermediate product range Storage tray and set list

The instrument set described as follows represents the basic equipment needed for implanting DMS Intermediate plates.

55-550-50-01
Storage tray



<i>DMS Intermediate Set</i>		
consisting of:		
Item No.	Quantity	Description
25-170-00-01	1	DMS Intermediate Set
Implants		
25-103-16-05 – 25-103-40-05	1	Cortical screw, 3.5 x 16-40 mm, self-tapping
Instruments		
22-368-24-07	1	Screwdriver, 2.5 mm, Allen
25-219-04-07	1	Depth gauge, small, 16 cm
25-210-25-07	1	Drill bit, 2.5 mm
25-228-00-05	1	DMS Intermediate guide wire
25-228-05-07	1	DMS Intermediate measuring sleeve
25-228-06-07	1	DMS Intermediate lag screw inserter
25-228-12-07	1	DMS Intermediate combo reamer
25-228-14-07	1	DMS Intermediate plate impactor
25-228-17-07	1	DMS Intermediate tap
25-228-20-07	1	DMS Intermediate lag screw connector
25-228-21-07	1	DMS Intermediate safety inserter for screwdriver
25-222-32-01	1	Drill bit guide for 3.5-mm screws
25-238-08-07	1	DMS T-handle
Storage tray		
55-550-50-01	1	DMS Intermediate storage tray for instruments and implants

Notice:

The DMS Intermediate plates, lag screws and – if needed – compression screws must be ordered separately for each patient.

References

- Dittel KK
Die Osteosynthese pertrochantärer Oberschenkelfrakturen mit einer neuen winkeladaptierbaren dynamischen Gleitlasche
57. Jahrestagung der Deutschen Gesellschaft für Unfallchirurgie e.V., Berlin
Abstract-Band **(1993)** 66
- Dittel KK
Dynamische Osteosynthese pertrochantärer Oberschenkel-frakturen mit einer winkeladaptierten Laschenplatte
Kongressband Osteosynthese International Erlangen 1993
Druckhaus Mayer Verlag Erlangen **(1994)** 95-102
- Dittel KK
An Innovative Method for Stabilization of Proximal Femur Fractures
Abstract-Book: Br. J. Surg. **(1994)** 81 Suppl 1, 126
- Dittel KK
Ein innovatives Stabilisierungsverfahren für proximale Femurfrakturen
Osteosynthese International (Kongressband),
Leuven University Press Verlag **(1995)** 193 - 198
- Dittel KK, Rapp M
Ein neues Prinzip zur Stabilisierung proximaler Femurfrakturen
Osteosynthese International **(1995)** 1, 46-54
- Dittel KK, Schier H, Rapp M
Actual results after fracture management at the proximal femur using a new dynamic angle-adaptable device
Osteosynthesis International; Oulu **(1995)**
Herausgeber: Jalovaara P., Vécsei V.; S. 319 - 324
- Molinar Min AM, Moselli M
La vite-placca DMS nell'osteosintesi delle fratture pertrocanteriche di femore
Minerva Ortop Traumatol **(1995)** 46, 329-330
- Rapp M, Dittel KK, Schier H
A new dynamic angle-adapted device – an innovative method for stabilizing proximal fractures of the femur
Br. J. Surg. **(1995)** 82 Suppl. 1, 119 (Abstract)
- Dittel KK, Schier H, Rapp M
Actual results after fracture management at the proximal femur using a new dynamic angle-adaptable device
Osteosynthese International (Kongressband) Oulu **(1996)** 319 - 324
- Dittel KK, Rapp M, Schier H
DMS: Un principio estabilizador innovador para fracturas proximales de fémur
Rev. S. and Traum. y Ort. **(1997)** 17-1, 73-81
- Dittel KK, Felenda MR
Operative Behandlung der Gelenk- und Schaftfrakturen
Thieme Verlag, Stuttgart – New York **(1998)**
- Dittel KK, Felenda MR
Spezielle Indikationsstellung zur Osteosynthese suprakondylärer Femurfrakturen mit der DMS in Wiederherstellungschirurgie des Kniegelenkes – Wandel in der Osteosynthesetechnik
Rahmanzadeh, Voigt, Trabhardt (Eichhorn Presseverlag) **(1998)** 148-153
- Dittel KK, Reinecke M
Erste klinische Erfahrungen mit einer zementfreien, extra metaphysär orientierten Schenkelhalsprothese
Hefte zu „Der Unfallchirurg“, Springer Verlag **(1998)** 718 - 719
- Rapp M, Dittel KK, Eberhard HJ, Miller WO
Eine innovative Methode zur Stabilisierung instabiler intertrochantärer Umkehrfrakturen des proximalen Femurs (31 A 3.3)
12. Internationaler Kongress „Osteosynthese International“ des Gerhard-Küntschers-Kreises e.V. Stuttgart
Abstract-Band 70 **(1998)**

- Rapp M, Miller WO, Dittel KK, Abendschein W
The variable angle compression hip system
12. Internationaler Kongress „Osteosynthese International“
des Gerhard-Küntschers-Kreises e.V., Stuttgart
Abstract-Band 156 **(1998)**
- Dittel KK, Felenda MR
Erste klinische Erfahrungen mit einer zementfreien
extra metaphysär orientierten Schenkelhalsprothese
Osteosynthese International **(1999)** Suppl. 2, 7, 164 - 167
- Rapp M, Eberhard HJ, Dittel KK, Miller WO
Eine innovative Methode zur Stabilisierung instabiler
intertrochantärer Umkehrfrakturen des proximalen
Femurs (31 A 3.3)
Jubiläumskongress „Osteosynthese International 1998“
des Gerhard-Küntschers-Kreises e.V., Stuttgart
Osteosynthese International **(1999)** 7 Suppl. 2, 47-51
- Ateschrang A
Die Prognose der medialen Schenkelhalsfraktur bei
kopferhaltender Osteosynthese
Inauguraldissertation Universität Tübingen **(2000)**
- Rapp M, Dittel KK, Felenda MR
Die Dynamische KLS-Martin-Schraube (D.M.S.) als
Implantatalternative zur Stabilisierung suprakondylärer
Femurfrakturen
35. Jahrestagung der Österreichischen Gesellschaft
für Unfallchirurgie, Salzburg/Österreich
Acta Chirurgica Austriaca **(2000)** 32, [Suppl. 161], 71 - 75
- Hajinpour MA
The Vari Angle Compression Hip System: A New Device
for the Treatment of Hip Fractures
Preliminary report of 114 cases
Journal of Trauma: **(2001)** 56 - 65
- Chaim SH et al
A Biomechanical Study of Femoral Neck Fracture Fixation
with the VHS Vari-Angle Hip Fixation System
Supplementum to the American Journal of Orthopaedics:
(2002) 22 - 24
- Dittel KK, Weise K
Komplikationsmanagement in der Traumatologie
Thieme Verlag, Stuttgart – New York **(2003)**
- Dittel KK, Rapp M
10 years of experience with the DMS (Dynamic KLS Martin
Screw) in the treatment of proximal and distal femoral
fractures.
Osteosynthese International Kongress **(2004)** Graz/Österreich,
Abstract-Band 41
- Dittel KK, Rapp M
Winkelstabile implants am proximalen Femur.
Akt Traumatol. **(2005)** 35, 155-162
- Ateschrang A, Dittel KK
The Dynamic KLS Martin Screw: an Alternative for
Intracapsular Femoral Neck Fractures?
Eur J Trauma Emerg Surg **(2007)** 33, 635-640
- Mishra Anil Kumra
Management of Intertrochanteric Fractures by using
Dynamic Hip Screw / Dynamic KLS Martin Screw
J. Orthopaedics **(2007)** 4 (2) e 40
- Dittel KK, Rapp M
The Double Dynamic KLS Martin Screw (DMS)
Steinkopferverlag Heidelberg, **(2008)**

KLS Martin Group

Karl Leibinger Medizintechnik GmbH & Co. KG
78570 Mühlheim · Germany
Tel. +49 7463 838-0
info@klsmartin.com

KLS Martin GmbH + Co. KG
79224 Umkirch · Germany
Tel. +49 7665 98 02-0
info@klsmartin.com

Stuckenbrock Medizintechnik GmbH
78532 Tuttlingen · Germany
Tel. +49 7461 165880
verwaltung@stuckenbrock.de

Rudolf Buck GmbH
78570 Mühlheim · Germany
Tel. +49 7463 99516-30
info@klsmartin.com

KLS Martin France SARL
68200 Mulhouse · France
Tel. +33 3 89 51 31 50
france@klsmartin.com

Martin Italia S.r.l.
20871 Vimercate (MB) · Italy
Tel. +39 039 605 67 31
italia@klsmartin.com

Martin Nederland/Marned B.V.
1271 AG Huizen · The Netherlands
Tel. +31 35 523 45 38
nederland@klsmartin.com

KLS Martin UK Ltd.
Reading RG1 3EU · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

Nippon Martin K.K.
Osaka 541-0046 · Japan
Tel. +81 6 62 28 90 75
nippon@klsmartin.com

KLS Martin L.P.
Jacksonville, FL 32246 · USA
Tel. +1 904 641 77 46
usa@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
121471 Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
201203 Shanghai · China
Tel. +86 21 2898 6611
china@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

**1923
2013
YEARS** | **SURGICAL
INNOVATION**

Gebrüder Martin GmbH & Co. KG
A company of the KLS Martin Group
KLS Martin Platz 1 · 78532 Tuttlingen · Germany
Postfach 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com

