


## IXOS® Radius Plating System

*Simply clever!*





***Our core competence** is hand surgery, a field where we can offer you much more than just standard treatment solutions for, say, distal radius fractures. Many of our products are intended to help you to achieve outstanding results in difficult, non-everyday situations as well. Products such as our ulnar head prosthesis (UHP) or the Flower Plate for mediocarpal partial arthrodesis are excellent examples of this.*

***Our objective** is to simplify hand surgery interventions through intelligent system solutions, helping you to achieve the best possible results in the interest of the patient. Working in close cooperation with well-known authors and their teams, we have translated new ideas into innovative products that are consistently being developed further in an ongoing process. The result is a wide range of high-quality systems that impress with their clever design along with easy and safe handling.*

*And what's more, we have never lost sight of the economic perspective and service needs of our customers.*

***We consider ourselves as a true partner** – to be relied upon for routine tasks and special challenges alike.*

## Table of contents

	Pages
IXOS® – product features	6-11
Indications and surgical techniques	12-13
■ Radius fracture	
Treatment with classical palmar plate	14-21
■ Radius fracture	
Treatment with palmar watershed-line plate	22-29
Product range	
■ IXOS® implants	30-35
■ IXOS® instruments	36-39
■ IXOS® storage system	40-41

IXOS®: *simply clever!*



## **IXOS®** *Radius Plating System* – simply clever!

The most frequent fractures encountered in surgical practice are radius fractures. Based on an extensive body of clinical experience and utilizing new technical possibilities afforded by the ever-advancing technical progress, multidirectional locking radius plates are frequently used nowadays for treating such conditions.

Based on the principle of multidirectional locking plate osteosynthesis, our goal was to treat nearly all types of distal radius fractures with an easy-to-use and clearly structured system.

In fact, IXOS® comes as a comprehensive and user-friendly radius plating system including palmar, dorsal as well as lateral plates. All plates are implanted with state-of-the-art smartDrive® screws. No more than four instruments are required for secure osteosynthesis.





## *Feature, Function* and Benefit



IXOS® radius plates are available in different designs to match proved treatment concepts. All plates are finished with the Dotize® surface coating. To facilitate identification, all palmar plates have been marked “P”, dorsal plates “D” and dorsolateral plates “DL”.

The latest generation of smartDrive® screws provides both standard and locking screws with double threads for the first time. In addition, all screws are equipped with atraumatic screw tip. Of course, the smartDrive® screws also exhibit the T8 with self-retaining function that has been established for decades. The product range is complemented by locking pins.

The screws/pins are color-coded to facilitate their application:


**Blue:** locking smartDrive® screws and pins

**Gold:** standard smartDrive® screws

## IXOS® – simply clever!


### Features

### Benefits



- Anatomical plate geometry
- Rounded atraumatic plate contour

- No need to bend plates
- Best possible embedding in soft tissue



- marLock locking
- Angulation within a range of  $\pm 15^\circ$
- Several times relockable
- Locking even without “heel piece”

- High degree of locking stability
- Best possible screw positioning
- Adjustment of screw position and easy metal removal
- Secure use of pins




Type II anodization

- 15% more fatigue resistance
- Smooth surface
- Risk of contact welding is minimized


- Slim plate design
- Delays adherence of tissue and bone ingrowth
- Easy metal removal

## smartDrive® screws




- Atraumatic screw tip
- Atraumatic screw head

- Secure bicortical anchoring with soft tissue preservation
- Maximum angulation without soft tissue irritations



- Double, self-tapping thread

- Reduces screw-in time by 50%



- T8 with self-retaining function

- Easy pick-up, insertion, tightening or removal of the screw

## *Feature, Function* and Benefit in Detail



IXOS® components are manufactured according to the latest findings. The 3D contour of the P4 and the P4 *Wave* can only be achieved by manufacturing them on state-of-the-art, computer-controlled 5-axle machines.

The following plate types are available:


- P2: This plate complies with the present industrial standard and complements the system with regard to economic aspects.
- P4: Based on the classic palmar treatment concept, the P4 exhibits unprecedented product features in this category for the first time.
- P4 *Wave*: A watershed-line plate of the latest generation.
- DL4: Anatomically pre-shaped plates for the dorsolateral treatment of radius fractures.
- D4: The system is complemented by anatomically designed plates for dorsal treatment.
- PU4: Additional ulnar plate for the treatment of distal ulnar neck and head fractures.



## A solution for every situation


### Features

### Benefits



- The appropriate plate for every radius fracture even in terms of economic aspects


- No second system is required during surgery
- The same instruments for all plates



P2

- The industrial standard

- Familiar technology at a reasonable price



P4

P4 *Wave*

- Watershed-line technology for the first time both in conventional and anatomical design


- Multidirectional locking but nevertheless prepositioned screws

- Integrated support for ideal screw positioning

- Best possible ulnar support in anatomical design

- Extra-long plates in conventional design

- For the treatment of complicated distal radius and shaft fractures



DL4

- Plates with small dimensions

- Allows dual-plate technology according to Rikli & Regazzoni

PU4

- Special ulnar plate

- For the treatment of distal ulnar neck and head fractures



D4

- Atraumatic frame plate
- A great number of multidirectional locking boreholes

- Dorsal support but nevertheless minimum soft tissue irritation
- High flexibility of treatment
- Easy fine adjustment with special bending pliers possible

## *Feature, Function* and Benefit




The KLS Martin Group is committed to developing color-coded instruments that can be handled easily and efficiently. The radius plating system comprises only 4 instruments. In order to comply with the specific requirements of the users, both the screwdriver and the depth gauge can be selected according to the specific personal preferences.

The storage concept already known from HBS2 has been adapted to the special requirements of radius treatment. Besides easy handling, the maintenance requirements were given top priority.

## IXOS® instruments and storage


### Features

### Benefits




- Color-coded instruments (purple)
  - smartDrive® screws Ø 2.5 mm
  - smartDrive® pins Ø 2.0 mm

- Easy identification of the respective instruments




- Single-part instruments with ergonomically shaped silicone handles

- Good tactile feedback
- No couplings that could lead to confusion
- No parts that could get lost




- Monoaxial drill guide
- Polyaxial drill guide

- Combined with prepositioned holes, allows short surgery periods
- Precise screw positioning in compliance with the maximum possible angulation of  $\pm 15^\circ$




- Screwdrivers and depth gauges are available in two different design variants

- Intuitive working with optimum ergonomics



- Stainless steel storage tray in honeycomb design combined with high performance plastic

- High degree of stability at low weight
- Good rinsing results due to large openings
- No water residues



- The instruments are arranged according to their sequence of use during the surgical procedure

- For easy and efficient instrumentation

## *Step by Step* optimal Fixation

### *Indications*

Acute distal radius fractures



Type A2  
Colles' fracture



Type B3  
Smith's fracture  
Reversed Barton's fracture



Type A3



Type C1




Type B1



Type C2



Type B2  
Barton's fracture



Type C3



## *Surgical Techniques*

### Radius fracture

*Treatment with classical palmar plate*

*Pages 14-21*



### Radius fracture

*Treatment with palmar watershed-line plate*

*Pages 22-29*






Source: Dr. Meyer, Saarbrücken

### *Patient positioning*

In addition to taking standard x-rays in A/P and sagittal planes in neutral position of the wrist, a high-resolution computer tomography should be conducted for the further clarification of intra-articular fractures.

If a central impression of the distal end of the radius is suspected, a carpal arthroscopy can additionally be conducted to clarify concomitant injuries and assess the reduction.




### *Patient positioning*

The patient is placed on the back. The hand that requires surgery is extended on the operating table in full supination of the forearm and under tourniquet control.


A cloth roll can be placed underneath the wrist as a reduction aid.





### *1. Henry's palmar approach*

A skin incision of 6 - 10 cm length is made on the distal forearm three centimeters proximal to the wrist. The flexor carpi radialis tendon (FCR) is exposed.




### *2. Exposure of the radius*

To obtain access to the pronator quadratus, the incision extends between the FCR and the radius artery.

The pronator quadratus is detached from the lateral edge of the radius to elevate an ulnar-based flap.


#### *Note:*

*If a post-traumatic sensitivity disorder in the median nerve distribution area or a latent carpal tunnel syndrome is detected, the incision is extended distally and the carpal tunnel is opened.*



### 3. *Exposure of the fracture*

The fragments and the fracture gap are exposed.



### 4. *Insertion of the plate*

The implant is selected according to the fracture pattern and the patient's anatomy.

The plate is placed centrally above the longitudinal axis in the direction of the distal edge of the radius.

The plate can be temporarily fixed with Ø 1.2 mm K-wires.


The first borehole is made into the slotted hole of the shaft using the monodirectional drill guide and the core hole drill (1 purple ring).



Core hole drill  
AO fitting Ø 2.0 mm



Drill guide  
monoaxial



### 5. Determination of the screw length

The correct screw length is determined using the depth gauge.



Source: Dr. Meyer, Saarbrücken

### 6. Insertion of the first shaft screw

The plate is fixed in the slotted hole with a golden standard screw.

The correct plate position and the anatomical reduction are checked under x-ray control in both planes.

It has to be ensured that the plate does not project over the watershed-line; this might cause irritation to the flexor tendons.

If necessary, the result has to be corrected and the plate displaced in longitudinal and/or lateral direction. The screw has to be loosened for this purpose.




Depth gauge  
AO principle



Depth gauge  
Single-hand principle




Screwdriver T8



#### 7. Insertion of another shaft screw

In order to be able to absorb optimally the forces in the shaft region during reduction, it is advisable to insert another shaft screw, preferably a blue locking screw, prior to the reduction, ensuring that the plate is positioned correctly.



#### 8. Fracture reduction


The tilted fracture is reduced under x-ray control. The bent hand is reduced by applying longitudinal traction combined with dorsal digital compression.

*Note:*

*If required, the fracture reduction can be fixed with K-wires.*



Screwdriver T8




#### 9a. Insertion of the distal screws

The first distal borehole is made using the polyaxial drill guide and the core hole drill (1 purple ring). The screw length is determined and a light blue locking screw is inserted.

##### Note:

The drill guide allows for a multidirectional angulation of  $\pm 15^\circ$ , so that fixed-angle locking is always ensured.

If the drill guide is not used, the permissible locking angle might be exceeded, which could lead to a lasting impairment of the angle stability.



#### 9b. Insertion of the distal screws

The monoaxial drill guide can be used as an alternative. This takes up the prepositioned angles in the plate.

##### Note:

When treating a fracture with the P2 plate, the polyaxial drill guide shall always be used for positioning the distal screws.




Core hole drill  
AO fitting Ø 2.0 mm

Drill guide  
polyaxial



Core hole drill  
AO fitting Ø 2.0 mm

Drill guide  
monoaxial



#### 10. Insertion of additional screws

All additional screws are inserted. For this purpose, drilling and measuring is performed as usual. The screws are positioned in the direction of the dorsal edge of the radius. If possible, the radially positioned screw should be inserted into the radial styloid process.


The subchondral position of the screws is checked under x-ray control.

If required, spongiosa or bone substitute can be inserted through the plate window.

#### Note:

*The screws in the first row should slightly be tilted proximally; by contrast, the screws in the second row should be tilted distally.*

*The subchondral screw arrangement according to the array principle provides optimum support for both the central region and the dorsal edge of the radius.*



#### 11. Wound closure

The wound is closed in layers.

Following the skin closure, a final x-ray image is taken.



Screwdriver T8





Source: Dr. Meyer, Saarbrücken

### *12. Postoperative treatment*

After completion of the surgery, a forearm splint is applied, which allows active finger movement.




Source: Prof. Liener, Stuttgart

### *Preoperative planning*

In addition to taking standard x-rays in A/P and sagittal planes in neutral position of the wrist, a high-resolution computer tomography should be conducted for the further clarification of intra-articular fractures.


If a central impression of the distal end of the radius is suspected, a carpal arthroscopy can additionally be conducted to clarify concomitant injuries and assess the reduction.



### *Patient positioning*


The patient is placed on the back. The hand that requires surgery is extended on the operating table in full supination of the forearm and under tourniquet control.

A cloth roll can be placed underneath the wrist as a reduction aid.



### *1. Palmar approach*

A skin incision of 6-10 cm length is made on the distal forearm three centimeters proximal to the wrist. The incision is extended distally at acute angle to the rascetta. The flexor carpi radialis tendon (FCR) is exposed.




### *2. Exposure of the radius*

To obtain access to the pronator quadratus, the approach extends between the FCR and the radial artery.

The pronator quadratus is detached from the lateral edge of the radius to elevate an ulnar-based flap.


#### *Note:*

*If a post-traumatic sensitivity disorder in the median nerve distribution area or a latent carpal tunnel syndrome is detected, the incision is extended distally and the carpal tunnel is opened.*



### 3. Exposure of the fracture

The fragments and the fracture gap are exposed.




### 4. Fracture reduction

The tilted fracture is reduced under x-ray control. The bent hand is reduced by applying longitudinal traction combined with dorsal digital compression.

*Note:*

*If required, the fracture reduction can be fixed with K-wires.*




### 5. Insertion of the plate

The implant is selected according to the fracture pattern and the patient's anatomy.

The plate is placed centrally above the longitudinal axis in the direction of the distal edge of the radius.

The plate can be temporarily fixed with Ø 1,2 mm K-wires. The K-wires can be positioned in such a way that the position of the plate to the distal radioulnar joint (DRUJ) as well as to the radiocarpal joint can be checked simultaneously.

The positioning of the plate will be controlled by image converter.



### 6. Insertion of the distal screws


The first borehole is made into the ulnar plate hole using the polyaxial drill guide and the core hole drill (1 purple ring).

*Note:*

*The drill guide allows for a multidirectional angulation of +/- 15°, so that fixed-angle locking is always ensured.*


*If the drill guide is not used, the permissible locking angle might be exceeded, which could lead to a lasting impairment of the angle stability.*





#### 7. Determination of the screw length

The correct screw length is determined using the depth gauge.



#### 8. Insertion of the distal screws

The monoaxial drill guide can be used as an alternative. This takes up the prepositioned angles in the plate.



Depth gauge  
AO principle



Depth gauge  
Single-hand principle




Core hole drill  
AO fitting Ø 2.0 mm



Drill guide  
monoaxial





### 9. Insertion of the screws

The plate is fixed with a golden standard screw.

All additional screws are inserted at a fixed angle. For this purpose, drilling and measuring is performed as usual. The screws are positioned in the direction of the dorsal edge of the radius. If possible, the radially positioned screw should be inserted into the radial styloid process.


The subchondral position of the screws is checked under x-ray control. The K-wires can subsequently be removed.

#### Note:

*The screws in the first row should be slightly tilted proximally; by contrast, the screws in the rear rows should be tilted distally. The subchondral screw arrangement according to the array principle provides optimum support for both the central region and the dorsal edge of the radius.*



Screwdriver T8



### 10. Insertion of the first shaft screw

The distal fragment is brought into the final position by pressing the proximal end of the plate in place.

The plate is fixed in the slotted hole with a standard screw. This allows for making fine adjustments to the distal fragment, if necessary.




Core hole drill  
AO fitting Ø 2.0 mm



Drill guide  
monoaxial




Screwdriver T8



*11. Insertion of the remaining shaft screws*

The remaining locking shaft screws are inserted. For this purpose, drilling and measuring is performed as usual.

If required, spongiosa or bone substitute can be inserted through the plate window.




*12. Wound closure*

The wound is closed in layers.  
Following the skin closure, a final x-ray is taken.



Screwdriver T8



Source: Prof. Liener, Stuttgart

### *13. Postoperative treatment*

After completion of the surgery, a forearm splint is applied, which allows active finger movement.

*Implants* **IXOS®**  
Palmar Radius Plates

P2


Length 43 mm  
Width 23 mm

Length 52 mm  
Width 27 mm

Length 71 mm  
Width 24.5 mm




26-912-10-09 right  
26-912-11-09 left

 = 2.0 mm




26-912-12-09 right  
26-912-13-09 left


 = 2.0 mm



26-912-14-09 right  
26-912-15-09 left

 = 2.0 mm

**STERILE**



#### Explanation of icons

- Titanium, Dotize®
- Packing unit
- Multidirectional locking
- Plate profile
- Implants in sterile packaging

## P4

Length 43 mm  
Width 23 mm

Length 95 mm  
Width 23 mm

*Drill guide block for*  
P4 plates with a width  
of 23 mm

Length 52 mm  
Width 27 mm

*Drill guide block for*  
P4 plates with a width  
of 27 mm



26-914-10-09 right  
26-914-11-09 left

= 2.0 mm



26-914-14-09 right  
26-914-15-09 left

= 2.2 mm

**STERILE**



26-950-50-07 right  
26-950-51-07 left



26-914-12-09 right  
26-914-13-09 left

= 2.0 mm



26-950-52-07 right  
26-950-53-07 left

## Implants **IXOS®**

### Dorsal and Dorsolateral Radius Plates as well as Ulnar Plate

#### *P4 Wave*

Length 43 mm  
Width 23 mm

Length 52 mm  
Width 23 mm

Length 70 mm  
Width 23 mm

Length 95 mm  
Width 23 mm

*Drill guide block for*  
P4 *Wave* plates with a width  
of 23 mm



1/1

26-914-20-09 right  
26-914-21-09 left

= 2.0 mm



1/1

26-914-22-09 right  
26-914-23-09 left

= 2.0 mm



1/1

26-914-26-09 right  
26-914-27-09 left

= 2.2 mm

**STERILE**



1/1

26-914-24-09 right  
26-914-25-09 left

= 2.2 mm

**STERILE**



1/1

26-950-54-07 right  
26-950-55-07 left





### IXOS® P4 Wave

- 0.5° ulnar, 8.5° distal
- 9.9° ulnar, 7.1° proximal
- 4.3° ulnar, 5.7° proximal
- 5.7° radial, 8.5° distal
- 1.4° ulnar, 1.4° distal
- 9.9° radial, 11.3° distal
- 9.9° radial, 10.7° distal
- 11.3° radial, 12.7° distal

### Explanation of icons

- Titanium, Dotize®
- Packing unit
- Multidirectional locking
- Plate profile

**STERILE** Implants in sterile packaging

## D4

Length 60 mm  
Width 32 mm

Length 60 mm  
Width 30 mm



26-914-30-09 right  
26-914-31-09 left

= 1.7 mm



26-914-33-09 right  
26-914-34-09 left

= 1.7 mm

## DL4

### Straight plate

Length 52 mm  
Width 7.5 mm



26-914-40-09

= 1.7 mm

### L plate

Length 43 mm  
Width 15 mm



26-914-42-09 right  
26-914-43-09 left

= 1.7 mm

## PU4 Ulnar plate

Length 53 mm  
Width 16 mm



26-914-41-09

= 1.7 mm

**STERILE**

*Implants* **IXOS®**  
Screws and Pins

*Screws*

Multidirectional  
locking screw

Ø 2.5 mm



Length	Art. no.
8 mm	26-905-08-09
9 mm	26-905-09-09
10 mm	26-905-10-09
11 mm	26-905-11-09
12 mm	26-905-12-09
13 mm	26-905-13-09
14 mm	26-905-14-09
15 mm	26-905-15-09
16 mm	26-905-16-09
17 mm	26-905-17-09
18 mm	26-905-18-09
19 mm	26-905-19-09
20 mm	26-905-20-09
22 mm	26-905-22-09
24 mm	26-905-24-09
26 mm	26-905-26-09
28 mm	26-905-28-09
30 mm	26-905-30-09


*Screws*

Standard  
cortical screw





Ø 2.5 mm



Length	Art. no.
8 mm	26-906-08-09
9 mm	26-906-09-09
10 mm	26-906-10-09
11 mm	26-906-11-09
12 mm	26-906-12-09
13 mm	26-906-13-09
14 mm	26-906-14-09
15 mm	26-906-15-09
16 mm	26-906-16-09
17 mm	26-906-17-09
18 mm	26-906-18-09
19 mm	26-906-19-09
20 mm	26-906-20-09
22 mm	26-906-22-09
24 mm	26-906-24-09
26 mm	26-906-26-09
28 mm	26-906-28-09
30 mm	26-906-30-09



#### Explanation of icons

-  Screw diameter 2.5 mm/  
Pin diameter 2.0 mm
-  Titanium, Dotize®
-  T-Drive
-  Packing unit

#### *Pins*

Multidirectional  
locking pin

Ø 2.0 mm



$\frac{1}{1}$



Length	Art. no.
14 mm	26-907-14-09
15 mm	26-907-15-09
16 mm	26-907-16-09
17 mm	26-907-17-09
18 mm	26-907-18-09
19 mm	26-907-19-09
20 mm	26-907-20-09
22 mm	26-907-22-09
24 mm	26-907-24-09
26 mm	26-907-26-09
28 mm	26-907-28-09
30 mm	26-907-30-09

## *Instruments* **IXOS®**

### *Standard instruments*



26-950-01-07  
*Drill guide  
polyaxial*  
15 cm / 6"



26-950-02-07  
*Drill guide  
monoaxial*  
15 cm / 6"



26-950-20-07  
*Core hole drill*  
AO fitting  
Ø 2.0 mm  
11 cm / 4 1/4"








26-950-06-07  
*Depth gauge*  
Single hand principle  
13 cm / 5"





#### Explanation of icons

-  Screw diameter 2.5 mm
-  Steel
-  Silicone
-  T-Drive
-  Packing unit

#### Standard instruments



1/2

26-950-18-07  
*Screwdriver T8*  
19 cm / 7 1/2"



1/2

26-950-13-07  
*K-wire dispenser*  
Ø 1.2 mm  
17.5 cm / 6 3/4"



1/2

22-627-12-05  
*K-wires*  
Ø 1.2 mm  
12 cm / 4 3/4"



*Instruments*  
**IXOS®**

*Optional instruments*



26-950-03-07  
*Drill guide*  
conventional  
15 cm / 6"



26-950-04-09  
*Joystick*  
cannulated  
Ø 2.0 mm  
41.5 mm



26-950-21-07  
*Core hole drill*  
AO fitting  
Ø 2.0 mm  
11 cm / 4 1/4"



STERILE



26-950-25-07  
*Gliding hole drill*  
Ø 2.5 mm  
11 cm / 4 1/4"








26-950-26-07  
*Gliding hole drill*  
Ø 2.5 mm  
11 cm / 4 1/4"




STERILE



#### Explanation of icons

-  Screw diameter 2.5 mm
-  Steel
-  Silicone
-  T-Drive
-  Packing unit

**STERILE**  Instrument in sterile packaging

### Optional instruments



1/2

26-950-05-07  
*Depth gauge*  
AO principle  
13 cm / 5"



1/2

26-950-19-07  
*Screwdriver T8*  
rotary  
19 cm / 7 1/2"



1/2

26-950-16-07  
*Screwdriver blade*  
T8/AO  
10 cm / 4"



1/2

26-950-37-07  
*Bending pliers*  
17.5 cm / 6 3/4"





## *Storage System* **IXOS®**

The IXOS® storage system consists of various modules.

All instruments that are absolutely imperative for a surgery are stored separately in the instrument tray.

Optional instruments such as gliding hole drills or bending pliers for the dorsal plates can also be stored separately in the storage cage. Furthermore, there is additional free storage space that can be used individually.

Depending on the version, the circular screw rack can accommodate 180 (single-sided) or 360 screws (double-sided), 5 pieces of each type and length. The double-sided screw rack additionally provides the opportunity to store locking pins. The circular screw rack can be stocked individually.

In addition to the standard inventory of instruments, the **IXOS® storage set no. 26-900-10-04** includes a selection of implants that are specifically tailored to the treatment of palmar radius fractures.



#### Storage system

55-910-33-04 Storage set consisting of:  
lid, instrument insert, storage cage, circular rack for plates, single-sided circular screw rack

55-910-34-04 Storage set consisting of:  
lid, instrument insert, storage cage, circular rack for plates, double-sided circular screw rack



55-910-59-04  
Lid

1  
unit(s)



55-910-38-04  
Instrument tray  
for storage

1  
unit(s)



55-910-36-04  
Storage cage

1  
unit(s)



55-910-35-04  
Circular rack  
for plates

1  
unit(s)



55-910-39-04  
smartDrive® Ø 2.5 mm  
circular screw rack  
for screws, single-sided

1  
unit(s)



55-910-37-04  
smartDrive® Ø 2.5 mm  
circular screw rack  
for screws, double-sided

1  
unit(s)

## IXOS® Radius Plating System – simply clever!



### *IXOS® clearly impressed the jury*

IXOS® symbolizes our new generation of hand and trauma surgery products and the new corporate design of this product line. This particularly applies to the newly designed instruments and the new storage concept. Therefore, we are absolutely delighted that IXOS® has won the **IF design award**.

The meaning of the term “design” is frequently but incorrectly reduced to the appearance of a product. In fact, the term has a much wider scope, including functional as well as aesthetic properties. Thus, “design” highlights features that give users exactly the added value they are looking for.



*BOS Driver*  
Battery operated screwdriver



*Angulus*  
Right angled screwdriver



*IXOS®*  
Radius plating system



*Limax®*  
Diode-pumped Nd:YAG laser



2005



2011



2012



2012

*Surgical Innovation is our passion.*

Among experts, the **IF design award** is considered the top international competition.

We have won the product design award now the fourth time with in the category medicine / health + care, but the first time with an implant system.

Altogether 1605 firms from more than 48 countries participated in the competition for this highly regarded award with 4322 products. The fact that this year only just about 30% of all applications were considered prize-worthy shows the rigor the 44-member jury applied to its decisions.

## KLS Martin Group

**Karl Leibinger Medizintechnik GmbH & Co. KG**  
78570 Mühlheim · Germany  
Tel. +49 7463 838-0  
info@klsmartin.com

**KLS Martin GmbH + Co. KG**  
79224 Umkirch · Germany  
Tel. +49 7665 98 02-0  
info@klsmartin.com

**Stuckenbrock Medizintechnik GmbH**  
78532 Tuttlingen · Germany  
Tel. +49 7461 165880  
verwaltung@stuckenbrock.de

**Rudolf Buck GmbH**  
78570 Mühlheim · Germany  
Tel. +49 7463 99516-30  
info@klsmartin.com

**KLS Martin France SARL**  
68200 Mulhouse · France  
Tel. +33 3 89 51 3150  
france@klsmartin.com

**Martin Italia S.r.l.**  
20864 Agrate Brianza (MB) · Italy  
Tel. +39 039 605 6731  
italia@klsmartin.com

**Martin Nederland/Marned B.V.**  
1271 AG Huizen · The Netherlands  
Tel. +31 35 523 45 38  
nederland@klsmartin.com

**KLS Martin UK Ltd.**  
Reading RG1 3EU · United Kingdom  
Tel. +44 1189 000 570  
uk@klsmartin.com

**Nippon Martin K.K.**  
Osaka 541-0046 · Japan  
Tel. +81 6 62 28 9075  
nippon@klsmartin.com

**KLS Martin L.P.**  
Jacksonville, FL 32246 · USA  
Tel. +1 904 641 7746  
usa@klsmartin.com

**KLS Martin do Brasil Ltda.**  
CEP 04.531-011 São Paulo · Brazil  
Tel.: +55 11 3554 2299  
brazil@klsmartin.com

**KLS Martin Australia Pty Limited**  
Artarmon NSW 2064 · Australia  
Tel.: +61 2 9439 5316  
australia@klsmartin.com

**KLS Martin Malaysia Sdn. Bhd.**  
10200 Penang · Malaysia  
Tel. +604 263 2566  
malaysia@klsmartin.com

**Gebrüder Martin GmbH & Co. KG**  
Representative Office  
121471 Moscow · Russia  
Tel. +7 499 792-76-19  
russia@klsmartin.com

**Gebrüder Martin GmbH & Co. KG**  
Representative Office  
201203 Shanghai · China  
Tel. +86 21 5820 6251  
china@klsmartin.com

**Gebrüder Martin GmbH & Co. KG**  
Representative Office  
Dubai · United Arab Emirates  
Tel. +971 4 454 16 55  
middleeast@klsmartin.com

**Gebrüder Martin GmbH & Co. KG**  
A company of the KLS Martin Group  
KLS Martin Platz 1 · 78532 Tuttlingen · Germany  
Postfach 60 · 78501 Tuttlingen · Germany  
Tel. +49 7461 706-0 · Fax +49 7461 706-193  
info@klsmartin.com · www.klsmartin.com