


OSTEOTEC

Silicone Finger Implant™

Manufactured from high tear resistant implant grade silicone



- Proven dorsal load transferring hinge design
- Eleven, evenly scaled sizes for comprehensive anatomical fit
- Simple, precise instrumentation
- Designed and manufactured by an independent British company



Suited for both MCP and
PIP joint replacement

Silicone Finger Implant™

Implant



Sizer



Instrument Set



Sizers and Sterilisation Tray



Implants

Catalogue No.	Size
OSTF-00	00
OSTF-0	0
OSTF-1	1
OSTF-2	2
OSTF-3	3
OSTF-4	4
OSTF-5	5
OSTF-6	6
OSTF-7	7
OSTF-8	8
OSTF-9	9

Instrument Sets

Catalogue No.	Description
OSTF-INS1	Set of eight Proximal and eight Distal Bone Rasps, Starter Awl and 11 Sizers
OSTF-INS2	Sizer Sterilisation Tray, Empty
OSTF-INS3	Instrument Tray, Empty

Sizers

Catalogue No.	Size
OSTF-S00	00
OSTF-S0	0
OSTF-S1	1
OSTF-S2	2
OSTF-S3	3
OSTF-S4	4
OSTF-S5	5
OSTF-S6	6
OSTF-S7	7
OSTF-S8	8
OSTF-S9	9
OSTF-SIZ1	Sizing Set, 1 of each, Blue, Non-Sterile



Individual Rasps also available on request

Caution: Federal law restricts this device to sale by or on the order of a physician



Osteotec Limited

9 Silver Business Park, Airfield Way, Christchurch, Dorset, BH23 3TA

Tel: +44 (0)1202 487 885 Fax: +44 (0)1202 487 886

Web: www.osteotec.co.uk Email: sales@osteotec.co.uk

