



MALLET FINGER

Bertrand Fritsch, MD
Centre de la Main - Toulon - France

METAHUS®



design: e-morceaux.fr / 33 6404 © AREX® 2015

www.arex.fr

HK² - METAHUS®

Single use External Fixator

One box = HK218 or HK215 or HK212

- 10 HK² units (same reference)

One HK² unit =

- 2 K-wires, same diameter (1.2 mm or 1.5 mm or 1.8 mm)

- 1 connector to join the 2 K-wires together

• Ref. **HK218**: DISTAL RADIUS - METACARPAL BONE

2 KW, Ø 1,8 mm

1 connector for Ø 1,8 mm



• Ref. **HK215**: CARPAL - METACARPAL - PHALANX

2 KW, Ø 1,5 mm

1 connector for Ø 1,5 mm



• Ref. **HK212**: CARPAL - PHALANX

2 KW, Ø 1,2 mm


1 connector for Ø 1,2 mm




1 MANOTTE® pliers 1.8LB



Bending



Cutting




Crimping

DISTAL RADIUS FRACTURE

Dr Gérard Hoel
University Hospital of St Denis - Reunion Island


HK²



5th METACARPAL BASE FRACTURE

Philippe Liverneaux, MD
University Hospital of Strasbourg - France


METAHUS®



BENNETH'S FRACTURE

Philippe Liverneaux, MD
University Hospital of Strasbourg - France

METAHUS®



5th METACARPAL NECK FRACTURE

Philippe Liverneaux, MD
University Hospital of Strasbourg - France

METAHUS®

