


LIGAMENTOTAXOR®


Adjustable distraction!

Quick and easy to use


2 possible assemblies: short or long blue plastic wires.



« Long » blue plastic wires



« Short » blue plastic wires



LIGAMENTOTAXOR® Kit description (double sealed sterile pack - single use)

- 2 « long » blue plastic wires
- 2 « short » blue plastic wires
- 2 KW Ø: 1.2 mm
- 1 KW Ø: 1.5 mm (« U » shape frame)
- 2 springs with open ends
- 4 KW tip protections
- 1 single use drill guide



DISTRIBUTED BY

CE 0499

Fabricant :
AREX®
BP20
F91125 Palaiseau cedex FRANCE

www.arex.fr

LIGAMENTOTAXOR®

ASSEMBLY



www.arex.fr/QR/LIGAMENTOTAXOR-2014.mp4

ADJUSTABLE DISTRACTION

